

DENSITY GRADIENT COLUMN APPARATUS QualiDG-1505

Density Gradient Column Apparatus QualiDG-1505

The QualiDG-1505 is a precision instrument **designed for highly accurate density determination** of solid plastics and resin-based materials. Ideal for quality control and material evaluation, it enables rapid detection of slight density differences in similar substances, including natural and synthetic fibers.

Using a linear gradient formed by two fully miscible liquids inside a glass column, the system suspends a calibrated glass float at the exact point where its density matches that of the surrounding solution—ensuring reliable and reproducible results.

Key Applications:

- Density analysis of high molecular weight compounds
- Testing natural and synthetic fibers
- Evaluation of resins and polymer materials
- Quality control through precise density comparison



FEATURES/ADVANTAGES

Precision Temperature Control

High-accuracy microcontroller with integrated heating/cooling maintains stable bath temperature within ±0.1°C, ensuring consistent results.

Low-Vibration Stirring System

External drive with internal spray stirring ensures uniform mixing while minimizing vibration for reliable gradient formation.

Durable, Transparent Bath Chamber

Made from high-clarity resin for excellent insulation, visibility, and vibration resistance.

Superior to Traditional Methods

Offers greater accuracy, ease of use, and a broader measurement range than immersion or float-sink techniques.

TECHNICAL SPECIFICATIONS

SPECIFICATIONS	VALUE
Model	QualiDG-1505
Measurement Range	0.80g to 2.70g/cm ³
Accuracy	0.0001 g/cm ³
Density Gradient Tube Dimensions	ID 50mm, length 900mm
Graduation Range	0-850mm, resolution: 1mm
Number of Test Stations	3 (or more)
Temperature Range	20°C to 30°C
Temperature Accuracy	±0.05° C
Standard Glass Floats	0.80g to 2.70g/cm³ accuracy 0.0001 g/cm³





Connect with us

Qualitest USA (Corporate Sales Office)

Toll-Free: 1.877.884.TEST (8378) | Fax: 954.697.8211 E-mail: info@qualitest-inc.com Corporate Address: 8201 Peters Rd., #1000, Plantation, FL 33324, USA.

Qualitest Singapore (ASIA PACIFIC Regional Office)

Tel: +65 6393 5480 | E-mail: singapore@qualitest-inc.com Address: 50 Raffles Place, Singapore Land Tower, Level 46, Singapore, 048623. Contact our **QualiTeam** today to find out how we can help your organization **select the most suitable testing solution** for your application, requirements, and budget.

Qualitest Canada & International

Tel: +1.905.944.9825 | Fax: +1.905.944.0304 E-mail: sales@qualitest-inc.com Address: 70 East Beaver Creek Rd., #9, Richmond Hill, Ontario L4B 3B2, Canada.

Qualitest Indonesia (Representative Office)

Tel: +62 21 2985 9522 | Fax: +62 21 2985 9889 E-mail: indonesia@qualitest-inc.com Address: One Pacific Place Level 11, Jl. Jend. Sudirman, Kav. 52-53, SCBD Area, Jakarta 12190, Indonesia.

Qualitest Latin America (Mexico and LATAM Region)

E-mail: ventas@qualitest-inc.com

Qualitest India

E-mail: india@qualitest-inc.com 15th Floor, Dev Corpora, Pokhran Road No.1, Eastern Express Highway, Thane, Maharashtra, Mumbai, 400601, India

Qualitest FZE (Regional GCC/ME Office)

Tel: +971 4 8819252 | Fax: +971 4 8819262 Email: gcc@qualitest-inc.com Address: Jafza One, BB 1610, Jebel Ali Free Zone, PO Box 261440, Dubai, UAE.





